

## Ni-NTA Agarose

Catalog: AS045

### ◆ Introduction

Ni-NTA Agarose is designed for purification of recombinant polyhistidine-tagged proteins from bacteria, insect, and mammalian cells. Ni-NTA Agarose exhibit high affinity and selectivity for polyhistidine-tagged fusion proteins, and allow protein purification under native or denaturing conditions.

### ◆ Characteristics

<b>Support</b>	<b>6% highly cross-linked agarose</b>
<b>Binding capacity</b>	<b>20mg His tag fusion protein/ml medium</b>
<b>Mean Bead size</b>	<b>100µm</b>
<b>Recommended flow rate</b>	<b>&lt;10ml/min</b>
<b>pH Limits</b>	<b>pH 3-12</b>
<b>working temperature</b>	<b>Room Temperature</b>
<b>Storage</b>	<b>4-8°C, 20% ethanol</b>

### ◆ Application

Purification of His-tagged fusion proteins.

### ◆ Note

- 1、 All the buffers should not contain EDTA, citrates or reductants such as  $\beta$ - mercaptoethanol.
- 2、 Sample can be dissolved by equilibration buffer (pH 5.5 ~ 8.5). The optimal pH for samples binding to the resin is 8.0-8.5.
- 3、 Samples should be centrifuged or passed through 0.45-µm filter to prevent clogging the column.
- 4、 Adding Tween-20 into equilibration buffer and elution buffer to a final concentration of 1% can reduce the non-specific binding.