


Validation Data



Flow cytometry: 1×10^6 293T (Transfection) cells were intracellularly-stained with ABflo® 594 Rabbit IgG isotype control (A23821,5 μ l/Test,left) or ABflo® 594 Rabbit anti-Human IL8 mAb (A24961,5 μ l/Test,right).

Flow cytometry: 1×10^6 THP-1 cells (untreated,left) and THP-1 cells (treated with 100ng/ml TPA for 24 hours and 5 μ g/mL LPS for 7 hours,right) were intracellularly-stained with ABflo® 594 Rabbit anti-Human IL8 mAb (A24961,5 μ l/Test).