

# Nestin Rabbit pAb

Catalog No.: A0484 **5 Publications**

## Basic Information

### Observed MW

280kDa

### Calculated MW

177kDa

### Category

Primary antibody

### Applications

WB,IF/ICC,FC

### Cross-Reactivity

Human

## Background

This gene encodes a member of the intermediate filament protein family and is expressed primarily in nerve cells.

## Recommended Dilutions

<b>WB</b>	1:500 - 1:1000
<b>IF/ICC</b>	1:20 - 1:50
<b>FC</b>	1:20 - 1:50

## Immunogen Information

### Gene ID

10763

### Swiss Prot

P48681

### Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

### Synonyms

Nbla00170; Nestin

## Contact

 | 400-999-6126

 | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification


### Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,pH7.3.

## Validation Data

---



Western blot analysis of various lysates using Nestin Rabbit pAb (A0484) at 1:1000 dilution.  
Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Basic Kit (RM00020).  
Exposure time: 5s.