

TRIB1 Rabbit pAb

Catalog No.: A10134 **1 Publications**

Basic Information

Observed MW

38kDa

Calculated MW

41kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Enables mitogen-activated protein kinase kinase binding activity and protein kinase inhibitor activity. Involved in several processes, including JNK cascade; negative regulation of lipopolysaccharide-mediated signaling pathway; and regulation of protein kinase activity. Located in cytoplasm and nucleus.

Recommended Dilutions

WB 1:1000 - 1:4000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

10221

Swiss Prot

Q96RU8

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

C8FW; GIG2; TRB1; GIG-2; SKIP1; TRB-1; TRIB1

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG


Purification

Affinity purification

Storage

Store at -20°C

Validation Data



Western blot analysis of various lysates using TRIB1 Rabbit pAb (A10134) at 1:1000 dilution.
Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 5s.