

METTL21C Rabbit pAb

Catalog No.: A10403 **1 Publications**

Basic Information

Observed MW

29kDa

Calculated MW

30kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Mouse, Rat

Background

Enables heat shock protein binding activity and protein-lysine N-methyltransferase activity. Involved in protein methylation. Located in nucleus. Part of protein-containing complex.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

196541

Swiss Prot

Q5VZV1

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

C13orf39; METTL21C

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG


Purification

Affinity purification

Storage

Store at -20°C

Validation Data



Western blot analysis of various lysates using METTL21C Rabbit pAb (A10403) at 1:1000 dilution.
Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 10s.