

RELL2 Rabbit pAb

Catalog No.: A10620

Basic Information

Observed MW

Refer to figures

Calculated MW

32kDa

Category

Primary antibody

Applications

ELISA

Cross-Reactivity

Human

Background

Predicted to enable collagen binding activity. Involved in positive regulation of p38MAPK cascade. Predicted to be located in extracellular matrix.

Recommended Dilutions

ELISA

Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

285613

Swiss Prot

Q8NC24

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

C5orf16; RELL2

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.