

B9D2 Rabbit pAb

Catalog No.: A10627

Basic Information

Observed MW

Refer to figures

Calculated MW

19kDa

Category

Primary antibody

Applications

IF/ICC,ELISA

Cross-Reactivity

Human, Mouse, Rat

Background

This gene encodes a B9 domain protein, which are exclusively found in ciliated organisms. The gene is upregulated during mucociliary differentiation, and the encoded protein localizes to basal bodies and cilia. Disrupting expression of this gene results in ciliogenesis defects.

Recommended Dilutions

IF/ICC 1:50 - 1:200

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

80776

Swiss Prot

Q9BPU9

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

MKS10; MKSR2; ICIS-1; JBTS34; MKSR-2; B9D2

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.