

SIGIRR Rabbit pAb

Catalog No.: A10665

Basic Information

Observed MW

Refer to figures

Calculated MW

46kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Predicted to enable NAD⁺ nucleosidase activity. Involved in negative regulation of DNA-binding transcription factor activity. Located in membrane.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

59307

Swiss Prot

Q6IA17

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

TIR8; IL-1R8; SIGIRR

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.