

UBN1 Rabbit pAb

Catalog No.: A10713

Basic Information

Observed MW

Refer to figures

Calculated MW

119kDa/122kDa/123kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Predicted to enable DNA binding activity. Predicted to be involved in DNA replication-independent chromatin assembly. Predicted to act upstream of or within keratinocyte differentiation. Located in bicellular tight junction. Is expressed in eye; limb; and limb mesenchyme. Orthologous to human UBN1 (ubinnuclein 1).

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

170644

Swiss Prot

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

1110029L11Rik; 2610108L02Rik; UBN1

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.