

CCAR1 Rabbit pAb

Catalog No.: A13595

Basic Information

Observed MW

132kDa

Calculated MW

133kDa

Category

Primary antibody

Applications

WB, IP, ELISA

Cross-Reactivity

Human, Mouse, Rat

Background

Enables RNA polymerase II cis-regulatory region sequence-specific DNA binding activity; nuclear receptor coactivator activity; and transcription corepressor activity. Involved in positive regulation of cell migration and positive regulation of cell population proliferation. Acts upstream of or within positive regulation of apoptotic process. Located in nuclear envelope lumen.

Recommended Dilutions

WB 1:500 - 1:1000

IP 0.5µg-4µg antibody for
200µg-400µg extracts of
whole cells

ELISA Recommended starting
concentration is 1 µg/mL.
Please optimize the
concentration based on your
specific assay requirements.

Immunogen Information

Gene ID

55749

Swiss Prot

Q8IX12

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

CCAR1

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG


Purification

Affinity purification


Storage

Store at -20°C

Validation Data



Immunoprecipitation analysis of 200 µg extracts of HeLa cells using 3 µg CCAR1 antibody (A13595). Western blot was performed from the immunoprecipitate using CCAR1 antibody (A13595) at a dilution of 1:500.



Western blot analysis of various lysates using CCAR1 Rabbit pAb (A13595) at 1:1000 dilution. Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (RM00020). Exposure time: 5s.