

KIF2B Rabbit pAb

Catalog No.: A13600

Basic Information

Observed MW

Refer to figures

Calculated MW

76kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Predicted to enable microtubule binding activity and microtubule motor activity. Involved in metaphase plate congression; microtubule depolymerization; and regulation of chromosome segregation. Located in intercellular bridge; mitotic spindle; and nucleolus.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

84643

Swiss Prot

Q8N4N8

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

KIF2B

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C