

FAIM2 Rabbit pAb

Catalog No.: A14110

Basic Information

Observed MW

Refer to figures

Calculated MW

35kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Mouse

Background

Involved in regulation of neuron apoptotic process. Acts upstream of or within negative regulation of extrinsic apoptotic signaling pathway via death domain receptors. Located in membrane raft.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

23017

Swiss Prot

Q9BWQ8

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

LFG; LFG2; NGP35; NMP35; TMBIM2; FAIM2

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.