

ZNF141 Rabbit pAb

Catalog No.: A14799

Basic Information

Observed MW

Calculated MW

55kDa

Category

Primary antibody

Applications

IHC-P, ELISA

Cross-Reactivity

Human, Mouse, Rat

Background

The protein encoded by this gene is a zinc finger protein that may be a tumor suppressor. Defects in this gene have been associated with autosomal recessive postaxial polydactyly type A.

Recommended Dilutions

IHC-P 1:50 - 1:200

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

7700

Swiss Prot

Q15928

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

D4S90; pHZ-44; ZNF141

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG


Purification

Affinity purification

Storage

Store at -20°C

Validation Data



Immunohistochemistry analysis of paraffin-embedded Mouse lung using ZNF141 Rabbit pAb (A14799) at dilution of 1:100 (40x lens). Microwave antigen retrieval performed with 0.01M PBS Buffer (pH 7.2) prior to IHC staining.