

KSR1 Rabbit pAb

Catalog No.: A15118

Basic Information

Observed MW

110kDa

Calculated MW

102kDa

Category

Primary antibody

Applications

WB,IF/ICC,ELISA

Cross-Reactivity

Human, Mouse, Rat

Background

Enables 14-3-3 protein binding activity; ATP binding activity; and protein C-terminus binding activity. Involved in positive regulation of MAPK cascade. Located in endoplasmic reticulum and membrane. Part of protein-containing complex. Implicated in breast adenocarcinoma. Biomarker of breast cancer.

Recommended Dilutions

WB 1:500 - 1:2000

IF/ICC 1:50 - 1:100

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

8844

Swiss Prot

Q8IVT5

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

KSR; RSU2; KSR1

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage


Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of various lysates using KSR1 Rabbit pAb (A15118) at 1:1000 dilution.
Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 10s.



Immunofluorescence analysis of C6 cells using KSR1 Rabbit pAb (A15118) at dilution of 1:100 (40x lens). Secondary antibody: Cy3-conjugated Goat anti-Rabbit IgG (H+L) (AS007) at 1:500 dilution. Blue: DAPI for nuclear staining.