

TCTN2 Rabbit pAb

Catalog No.: A15191

Basic Information

Observed MW

Calculated MW

77kDa

Category

Primary antibody

Applications

ELISA

Cross-Reactivity

Human

Background

This gene encodes a type I membrane protein that belongs to the tectonic family. Studies in mice suggest that this protein may be involved in hedgehog signaling, and essential for ciliogenesis. Mutations in this gene are associated with Meckel syndrome type 8. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

Recommended Dilutions

ELISA

Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

79867

Swiss Prot

Q96GX1

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

MKS8; TECT2; JBTS24; C12orf38; TCTN2

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.