

# LSM11 Rabbit pAb

Catalog No.: A15210

## Basic Information

### Observed MW

### Calculated MW

40kDa

### Category

Primary antibody

### Applications

IF/ICC,ELISA

### Cross-Reactivity

Human

## Background

Enables U7 snRNA binding activity. Involved in positive regulation of G1/S transition of mitotic cell cycle. Located in nuclear body. Part of U7 snRNP and telomerase holoenzyme complex. Implicated in Aicardi-Goutieres syndrome.

## Recommended Dilutions

**IF/ICC** 1:50 - 1:200

**ELISA** Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

## Immunogen Information

### Gene ID

134353

### Swiss Prot

P83369

### Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

### Synonyms

AGS8; LSM11

## Contact

☎ | 400-999-6126

✉ | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

🌐 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.