

XPO4 Rabbit pAb

Catalog No.: A15503

Basic Information

Observed MW

Refer to figures

Calculated MW

130kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

XPO4 belongs to a large family of karyopherins (see MIM 602738) that mediate the transport of proteins and other cargo between the nuclear and cytoplasmic compartments (Lipowsky et al., 2000 [PubMed 10944119]).

Recommended Dilutions

WB 1:200 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

64328

Swiss Prot

Q9C0E2

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

exp4; XPO4

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.