

ZCCHC10 Rabbit pAb

Catalog No.: A16539

Basic Information

Observed MW

Calculated MW

21kDa

Category

Primary antibody

Applications

ELISA

Cross-Reactivity

Mouse

Background

Predicted to enable nucleic acid binding activity and zinc ion binding activity.

Recommended Dilutions

ELISA

Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

54819

Swiss Prot

Q8TBK6

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

ZCCHC10

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.