

ANGPTL7 Rabbit pAb

Catalog No.: A17606

Basic Information

Observed MW

Refer to figures

Calculated MW

40kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Enables identical protein binding activity. Involved in negative regulation of vasculature development involved in avascular cornea development in camera-type eye and regulation of extracellular matrix organization. Located in extracellular region.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

10218

Swiss Prot

O43827

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

AngX; CDT6; dj647M16.1; ANGPTL7

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.