

SLC25A2 Rabbit pAb

Catalog No.: A17788

Basic Information

Observed MW

Refer to figures

Calculated MW

33kDa

Category

Primary antibody

Applications

IF/ICC,ELISA

Cross-Reactivity

Mouse

Background

This intronless gene encodes a protein that localizes to the mitochondrial inner membrane and likely functions as a transporter of small molecules such as ornithine. This gene is located between the protocadherin beta and gamma gene clusters on chromosome 5.

Recommended Dilutions

IF/ICC 1:50 - 1:200

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

83884

Swiss Prot

Q9BXI2

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

ORC2; ORNT2; SLC25A2

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.