

EIL1 Rabbit pAb

Catalog No.: A18497

Basic Information

Observed MW

Refer to figures

Calculated MW

66kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Arabidopsis thaliana

Background

ethylene-insensitive3-like1 (EIL1)

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

817247

Swiss Prot

Q9SLH0

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

AtEIL1; ETHYLENE-INSENSITIVE3-like 1; T20P8.10; T20P8_10; EIL1

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.