

IAA3 Rabbit pAb

Catalog No.: A18859

Basic Information

Observed MW

Refer to figures

Calculated MW

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Arabidopsis thaliana

Background

SHY2/IAA3 regulates multiple auxin responses in roots. It is induced rapidly by IAA, and has been shown to be phosphorylated by oat phytochrome A in vitro.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

839570

Swiss Prot

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

F19P19.32; F19P19_32; IAA3; indole-3-acetic acid inducible 3; SHORT HYPOCOTYL 2

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.