

# MSI4/FVE Rabbit pAb

Catalog No.: A20204

## Basic Information

### Observed MW

Refer to figures

### Calculated MW

113kDa

### Category

Primary antibody

### Applications

WB, ELISA

### Cross-Reactivity

Arabidopsis thaliana

## Background

Controls flowering.

## Recommended Dilutions

**WB** 1:1000 - 1:2000

**ELISA** Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

## Immunogen Information

### Gene ID

816471

### Swiss Prot

P09874

### Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

### Synonyms

ACG1; ATMSI4; F3P11.12; F3P11\_12; MSI4; MULTICOPY SUPPRESSOR OF IRA1 4; NFC04; NFC4; NUCLEOSOME/CHROMATIN ASSEMBLY FACTOR GROUP C; WD-40 REPEAT PROTEIN MSI4; MSI4/FVE

## Contact

 | 400-999-6126

 | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.