

SOS2 Rabbit pAb

Catalog No.: A20292

Basic Information

Observed MW

Refer to figures

Calculated MW

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Arabidopsis thaliana

Background

encodes a member of the CBL-interacting protein kinase family, is a regulatory component controlling plant potassium nutrition

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

833502

Swiss Prot

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

ATSOS2; CBL-INTERACTING PROTEIN KINASE 24; CIPK24; K21B8.3; K21B8_3; SALT OVERLY SENSITIVE 2; SNF1-RELATED PROTEIN KINASE 3.11; SNRK3.11; SOS2

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.