

WDR72 Rabbit pAb

Catalog No.: A20392 **1 Publications**

Basic Information

Observed MW

Refer to figures

Calculated MW

123kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse, Rat

Background

This gene encodes a protein with eight WD-40 repeats. Mutations in this gene have been associated with amelogenesis imperfecta hypomaturation type 2A3. Alternative splicing results in multiple transcript variants.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

256764

Swiss Prot

Q3MJ13

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

AI2A3; WDR72

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C