

FER Rabbit pAb

Catalog No.: A20542

Basic Information

Observed MW

Refer to figures

Calculated MW

98kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Arabidopsis thaliana

Background

Encodes a synergid-expressed, plasma-membrane localized receptor-like kinase that accumulates asymmetrically in the synergid membrane at the filiform apparatus and mediates male-female gametophyte interactions during pollen tube reception. Also involved in powdery mildew infection.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

824318

Swiss Prot

Q9SCZ4

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

FERONIA; FER

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.