# OsISA1 Rabbit pAb

Catalog No.: A20853



### **Basic Information**

### **Observed MW**

Refer to figures

### **Calculated MW**

### Category

Primary antibody

## **Applications**

WB,ELISA

### **Cross-Reactivity**

Oryza sativa

## **Background**

## **Recommended Dilutions**

**WB** 1:500 - 1:2000

**ELISA** 

Recommended starting concentration is 1 µg/mL.
Please optimize the concentration based on your specific assay

requirements.

## **Immunogen Information**

Gene ID Swiss Prot

4346068

#### **Immunogen**

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

### **Synonyms**

ISA; SU1; ISA1; OsISA1; isoamylase

### **Contact**

6		400-999-6126
$\bowtie$		cn.market@abclonal.com.cn
<u></u>	Π	www.abclonal.com.cn

### **Product Information**

SourceIsotypePurificationRabbitIgGAffinity purification

### **Storage**

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.