

PE Mouse IgG1 isotype control

Catalog No.: A22878 **1 Publications**

Basic Information

Observed MW

Refer to figures

Calculated MW

Category

Primary antibody

Applications

FC

Cross-Reactivity

CloneNo number

AMC0553

Conjugate

PE. Ex:565nm. Em:574nm.

Background

Recommended Dilutions

FC 10 μ L/Test

Immunogen Information

Gene ID

Swiss Prot

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Mouse

Isotype

Mouse IgG1

Purification


Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.1% Proclin300,0.2% BSA,pH7.2.

Validation Data



Flow cytometry: 1×10^6 Human PBMC were surface-stained with PE Mouse anti-Human CD1a mAb (A22793, orange line) or PE Mouse IgG1 isotype control mAb (A22878, blue line).