# LOC\_Os07g41310 Rabbit pAb

Catalog No.: A23143



## **Basic Information**

#### **Observed MW**

Refer to figures

### **Calculated MW**

50kDa

### Category

Primary antibody

### **Applications**

WB,ELISA

#### **Cross-Reactivity**

Oryza sativa

# **Background**

# **Recommended Dilutions**

**WB** 1:500 - 1:2000

**ELISA** 

Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

# **Immunogen Information**

**Gene ID Swiss Prot** 4343844 Q6Z4G8

#### **Immunogen**

Synthetic peptide. This information is considered to be commercially sensitive.

# **Synonyms**

BC1L6; OsJ\_25037; LOC\_Os07g41310

### **Contact**

<u>a</u>		400-999-6126
$\bowtie$		cn.market@abclonal.com.cn
$\overline{a}$	ī	www.ahclonal.com.cn

### **Product Information**

**Purification Source** Isotype Rabbit Affinity purification IgG

#### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.