

BRRF2 Rabbit mAb

Catalog No.: A23835 **Recombinant**

Basic Information

Observed MW

72kDa

Calculated MW

56kDa

Category

Primary antibody

Applications

WB,IF/ICC,ELISA

Cross-Reactivity

Epstein-Barr virus

CloneNo number

ARC61624

Background

BRRF2 is a capsid protein, but it can be secreted into the bloodstream, giving it the opportunity to be a diagnostic or immunotherapy target. The Epstein-Barr virus BRRF2 gene product is involved in viral progeny production.

Recommended Dilutions

WB 1:2500 - 1:10000

IF/ICC 1:100 - 1:400

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

3783729

Swiss Prot

P03210

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

Tegument protein BRRF2; BRRF2

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol and 0.05% BSA, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.

Validation Data

