

Smad9 Rabbit pAb

Catalog No.: A24036

Basic Information

Observed MW

Calculated MW
52kDa

Category
Primary antibody

Applications
WB

Cross-Reactivity
Human

Background

The protein encoded by this gene is a member of the SMAD family, which transduces signals from TGF-beta family members. The encoded protein is activated by bone morphogenetic proteins and interacts with SMAD4. Two transcript variants encoding different isoforms have been found for this gene.

Recommended Dilutions

Immunogen Information

Gene ID
4093

Swiss Prot
O15198

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

PPH2; MADH6; MADH9; SMAD8; SMAD8A; SMAD8B; SMAD8/9; Smad9

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source
Rabbit

Isotype
IgG

Purification
Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.