

# Cytokeratin(HMW) Mouse mAb

Catalog No.: A24141

## Basic Information

---

**Observed MW**

**Calculated MW**

**Category**

Primary antibody

**Applications**

IHC-P

**Cross-Reactivity**

Human

**CloneNo number**

AMC0667

## Background

---

## Recommended Dilutions

---

IHC-P 1:100 - 1:200

## Immunogen Information

---

**Gene ID**

**Swiss Prot**

**Immunogen**

This information is considered to be commercially sensitive.

**Synonyms**

## Contact

---

 | 400-999-6126

 | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

---

**Source**

Mouse

**Isotype**

IgG

**Purification**

Affinity purification


**Storage**

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol and 0.05% BSA, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.

## Validation Data

---



Immunohistochemistry analysis of paraffin-embedded Human prostate tissue using Cytokeratin(HMW) Mouse mAb (A24141) at a dilution of 1:100 (20x lens). High pressure antigen retrieval performed with 0.01M Tris-EDTA Buffer (pH 9.0) prior to IHC staining.