

ABflo® 488 Rabbit anti-Human CD72 mAb

Catalog No.: A24236

Basic Information

Observed MW

Calculated MW

40kDa

Category

Primary antibody

Applications

FC

Cross-Reactivity

Human

CloneNo number

ARC62240

Conjugate

ABflo® 488. Ex:491nm. Em:516nm.

Background

Predicted to enable signaling receptor binding activity. Predicted to be involved in cell adhesion. Is integral component of plasma membrane.

Recommended Dilutions

FC 5 µl per 10⁶ cells in 100 µl volume

Immunogen Information

Gene ID

971

Swiss Prot

P21854

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

CD72; CD72b; LYB2; CD72 molecule

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at 2-8°C. Avoid freeze.

Buffer: PBS containing 0.2% BSA, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.

Validation Data



Flow cytometry: 1X10⁶ A549 cells(negative control,left) and Daudi cells(right) were surface-stained with ABflo® 488 Rabbit anti-Human CD72 mAb(A24236,5 µl/Test,orange line) or ABflo® 488 Rabbit IgG isotype control (A22069,5 µl/Test,blue line).Non-fluorescently stained cells were used as blank control (red line).

Flow cytometry:1X10⁶ Daudi cells were surface-stained with ABflo® 488 Rabbit IgG isotype control (A22069,5 µl/Test,left) or ABflo® 488 Rabbit anti-Human CD72 mAb(A24236,5 µl/Test,right).