

LRP4 Rabbit pAb

Catalog No.: A25225

Basic Information

Observed MW

Refer to figures

Calculated MW

212kDa

Category

Primary antibody

Applications

IF-P, ELISA

Cross-Reactivity

Human

Background

This gene encodes a member of the low-density lipoprotein receptor-related protein family. The encoded protein may be a regulator of Wnt signaling. Mutations in this gene are associated with Cenani-Lenz syndrome.

Recommended Dilutions

IF-P 1:50 - 1:200

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

4038

Swiss Prot

O75096

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

CLSS; CMS17; LRP-4; LRP10; MEGF7; SOST2; LRP4

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG


Purification

Affinity purification

Storage

Store at -20°C

Validation Data



Immunofluorescence analysis of paraffin-embedded Human skeletal muscle tissue using LRP4 Rabbit pAb(A25225) at a dilution of 1:100 (40x lens). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L)(AS007) at 1:500 dilution. Blue: DAPI for nuclear staining. Perform high pressure antigen retrieval with 0.01 M citrate buffer (pH 6.0) prior to IF staining.