

ABflo® 594 Mouse IgG 2a isotype control

Catalog No.: A25263

Basic Information

Observed MW

Calculated MW

Category

Primary antibody

Applications

FC

Cross-Reactivity

CloneNo number

AMC50006

Conjugate

ABflo® 594. Ex:588nm. Em:604nm.

Background

Recommended Dilutions

FC 5 µl per 10⁶ cells in
100 µl volume

Immunogen Information

Gene ID

Swiss Prot

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Mouse

Isotype

IgG

Purification


Affinity purification

Storage

Store at 2-8°C. Avoid freeze.

Buffer: PBS with 0.03% proclin300,pH7.3.

Validation Data



Flow cytometry: 1×10^6 Raji cells (negative control, left) and Human PBMC (right) were surface-stained with ABflo® 594 Mouse anti-Human CD3 mAb (A25188, 5 μ l/Test, orange line) or ABflo® 594 Mouse IgG 2a isotype control (A25263, 5 μ l/Test, blue line). Non-fluorescently stained cells were used as blank control (red line).