

ZNF217 Rabbit pAb

Catalog No.: A27151

Basic Information

Observed MW

Calculated MW

115kDa

Category

Primary antibody

Applications

WB,ELISA

Cross-Reactivity

Human

Background

Enables DNA-binding transcription repressor activity, RNA polymerase II-specific and RNA polymerase II cis-regulatory region sequence-specific DNA binding activity. Involved in negative regulation of transcription by RNA polymerase II. Located in mitochondrion and nuclear speck. Part of histone deacetylase complex.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

7764

Swiss Prot

O75362

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

ZABC1

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.