

PE Rabbit anti-Mouse CD6

Catalog No.: A27330

Basic Information

Observed MW

Calculated MW

64kDa/68kDa/72kDa

Category

Primary antibody

Applications

FC

Cross-Reactivity

Mouse

CloneNo number

ARC66616

Conjugate

PE. Ex:565nm. Em:574nm.

Background

Predicted to enable several functions, including identical protein binding activity; lipopolysaccharide binding activity; and lipoteichoic acid binding activity. Predicted to be involved in several processes, including heterophilic cell-cell adhesion via plasma membrane cell adhesion molecules; lipopolysaccharide-mediated signaling pathway; and positive regulation of macromolecule metabolic process. Predicted to act upstream of or within cell adhesion. Predicted to be located in immunological synapse. Predicted to be spanning component of plasma membrane. Predicted to be part of T cell receptor complex. Predicted to be active in external side of plasma membrane. Orthologous to human CD6 (CD6 molecule).

Recommended Dilutions

FC ≤0.125 µg per million cells
in 100 µl volume

Immunogen Information

Gene ID

12511

Swiss Prot

Q61003

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG


Purification

Affinity purification

Storage

Store at -20°C

Validation Data



Flow cytometry: 1×10^6 NIH/3T3 cells (negative control, left) and C57BL/6 mouse splenocytes (right) were surface-stained with PE Rabbit anti-Mouse CD6 mAb (A27330, 0.125 μ g, orange line) or PE Rabbit IgG isotype control (A24172, 5 μ l/Test, blue line). Non-fluorescently stained cells were used as blank control (red line).

Flow cytometry: 1×10^6 C57BL/6 mouse splenocytes were surface-stained with PE Rabbit IgG isotype control (A24172, 5 μ l/Test, left) or PE Rabbit anti-Mouse CD6 mAb (A27330, 0.125 μ g, right).