

PE/Cyanine7 anti Human CD223/LAG-3 mAb

Catalog No.: A27492

Basic Information

Observed MW

Calculated MW

57kDa

Category

Primary antibody

Applications

FC

Cross-Reactivity

Human

CloneNo number

AMC0821

Conjugate

PE-Cy7. Ex:565nm. Em:778nm.

Background

Lymphocyte-activation protein 3 belongs to Ig superfamily and contains 4 extracellular Ig-like domains. The LAG3 gene contains 8 exons. The sequence data, exon/intron organization, and chromosomal localization all indicate a close relationship of LAG3 to CD4.

Recommended Dilutions

FC 5 μ l per 10^6 cells in 100 μ l volume

Immunogen Information

Gene ID

3902

Swiss Prot

P18627

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

CD223

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Mouse

Isotype

IgG1, kappa

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze.

Buffer: PBS with 0.09% Sodium azide, 0.2% BSA, pH7.2.

Validation Data

