

PE Rat anti Mouse CD206/MMR mAb

Catalog No.: A27501 **1 Publications**

Basic Information

Observed MW

Calculated MW

165kDa

Category

Primary antibody

Applications

FC (intra)

Cross-Reactivity

Mouse

CloneNo number

ATC0035

Conjugate

PE. Ex:565nm. Em:574nm.

Background

Enables mannose binding activity and transmembrane signaling receptor activity. Acts upstream of or within cellular response to interferon-gamma; cellular response to interleukin-4; and cellular response to lipopolysaccharide. Located in cell surface. Is expressed in several structures, including hemolymphoid system; intestine; liver; lung; and yolk sac. Orthologous to human MRC1 (mannose receptor C-type 1).

Recommended Dilutions

FC 0.125 µg per 10⁶ cells in
100 µl volume

Immunogen Information

Gene ID

17533

Swiss Prot

Q61830

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

MR; CD206

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rat

Isotype

Rat / IgG2b, kappa

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze.

Buffer: PBS with 0.09% Sodium azide, pH7.2.

Validation Data

