

ABflo® 450 Rat anti-Mouse CD206/MMR mAb

Catalog No.: A27626

Basic Information

Observed MW

Calculated MW

165kDa

Category

Primary antibody

Applications

FC

Cross-Reactivity

Mouse

CloneNo number

ATC0058

Conjugate

ABflo® 450. Ex:406nm. Em:445nm.

Background

Enables D-mannose binding activity and transmembrane signaling receptor activity. Acts upstream of or within cellular response to interleukin-4; cellular response to lipopolysaccharide; and cellular response to type II interferon. Located in cell surface. Is expressed in several structures, including hemolymphoid system; intestine; liver; lung; and yolk sac. Orthologous to human MRC1 (mannose receptor C-type 1).

Recommended Dilutions

FC ≤0.5 µg per million cells
in 100 µl volume

Immunogen Information

Gene ID

17533

Swiss Prot

Q61830

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

MR; CD206

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rat

Isotype

Rat / IgG2b, kappa

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze.

Buffer: PBS with 0.09% Sodium azide, pH7.2.