

# ABflo® 594 Rat anti-Mouse Ly-6G/Ly-6C (Gr-1)mAb

Catalog No.: A27642

## Basic Information

### Observed MW

### Calculated MW

14kDa

### Category

Primary antibody

### Applications

FC

### Cross-Reactivity

Mouse

### CloneNo number

ATC0033

### Conjugate

ABflo® 594. Ex:588nm. Em:604nm.

## Background

Located in external side of plasma membrane.

## Recommended Dilutions

**FC** ≤0.125 µg per million cells in 100 µl volume

## Immunogen Information

### Gene ID

546644/17067

### Swiss Prot

P35461/P0CW02

### Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

### Synonyms

Gr1; Gr-1; Ly-6G

## Contact

☎ | 400-999-6126

✉ | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

🌐 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Rat

### Isotype

IgG

### Purification

Affinity purification


### Storage

Store at 2-8°C. Avoid freeze.

Buffer: PBS with 0.09% Sodium azide, 0.2% BSA, pH7.3.

## Validation Data

---



Flow cytometry:  $1 \times 10^6$  C57BL/6 bone marrow cells were surface-stained with ABflo® 594 Rat IgG isotype control (5 µl/Test, left) or ABflo® 594 Rat anti-Mouse Ly-6G/Ly-6C (Gr-1) mAb (A27642, 0.125 µg, right).