

ABflo® 594 Rat anti-Mouse Ly-6G mAb

Catalog No.: A27647

Basic Information

Observed MW

Calculated MW

Category

Primary antibody

Applications

FC

Cross-Reactivity

Mouse

CloneNo number

ATC0030

Conjugate

ABflo® 594. Ex:588nm. Em:604nm.

Background

Located in external side of plasma membrane.

Recommended Dilutions

FC ≤0.25 µg per million cells
in 100 µl volume

Immunogen Information

Gene ID

546644

Swiss Prot

P35461

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

Gr1; Gr-1; Ly-6G

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rat

Isotype

Rat / IgG2a, kappa

Purification


Affinity purification

Storage

Store at 2-8°C. Avoid freeze.

Buffer: PBS with 0.09% Sodium azide, 0.2% BSA, pH7.3.

Validation Data



Flow cytometry: 1X10⁶ C57BL/6 mouse bone marrow cells were surface-stained with ABflo® 488 Rabbit anti-Mouse CD11b mAb (A24094, 5 µl/Test) and ABflo® 594 Rat IgG isotype control (5 µl/Test, left) or ABflo® 594 Rat anti-Mouse Ly6G mAb (A27647, 5 µl/Test, right).