

PE anti-Human Bcl-2 mAb

Catalog No.: A27722

Basic Information

Observed MW

Calculated MW
26kDa

Category
Primary antibody

Applications
FC

Cross-Reactivity
Human

CloneNo number
AMC0840

Conjugate
PE. Ex:565nm. Em:574nm.

Background

This gene encodes an integral outer mitochondrial membrane protein that blocks the apoptotic death of some cells such as lymphocytes. Constitutive expression of BCL2, such as in the case of translocation of BCL2 to Ig heavy chain locus, is thought to be the cause of follicular lymphoma. Alternative splicing results in multiple transcript variants.

Recommended Dilutions

FC 5 μ l per 10^6 cells in 100 μ l
volume

Immunogen Information

Gene ID	Swiss Prot
596	P10415

Immunogen
Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms
Bcl-2; PPP1R50

Contact

	400-999-6126
	cn.market@abclonal.com.cn
	www.abclonal.com.cn

Product Information

Source	Isotype	Purification
Mouse	Mouse IgG1	Affinity purification

Storage
Store at 4°C. Avoid freeze.
Buffer: PBS with 0.09% Sodium azide, 0.2% BSA, pH7.2.