

Magnetic Beads-conjugated Rabbit anti-Human TMEM192 mAb [MJFF-G-20C1]

Catalog No.: A27784

Basic Information

Observed MW

31 kDa

Calculated MW

31 kDa

Category

Primary antibody

Applications

IP

Cross-Reactivity

Human

CloneNo number

ARC5176-Magnetic Beads

Conjugate

Magnetic Beads

Recommended Dilutions

IP 0.5 µg - 4 µg antibody for
400µg-600µg extracts of
whole cells

Background

Enables protein homodimerization activity. Located in several cellular components, including late endosome; lysosomal membrane; and perinuclear region of cytoplasm.

Immunogen Information

Gene ID

201931

Swiss Prot

Q8IY95

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

TMEM192

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG


Purification

Affinity purification

Storage

Store at 2-8°C. Avoid freeze / thaw cycles.
Buffer: PBS with 0.09% Sodium azide, pH7.3.

Validation Data



Immunoprecipitation of Human TMEM192 from 600 µg extracts of U-251 MG cells was performed using 1 µg of Magnetic Beads-conjugated Rabbit anti-Human TMEM192 mAb [MJFF-G-20C1](A27784). Magnetic Beads-conjugated Rabbit IgG isotype control (AC047) was used to precipitate the Control IgG sample. IP samples were eluted with 1x Laemmli Buffer. The Input lane represents 10% of the total input. Western blot analysis of immunoprecipitates was conducted using TMEM192 Rabbit mAb (A21165) at a dilution of 1:1000.