

PE Mouse anti-Human TCR V62 mAb

Catalog No.: A27981

Basic Information

Observed MW

Calculated MW

19kDa

Category

Primary antibody

Applications

FC

Cross-Reactivity

Human

CloneNo number

AMC10019

Conjugate

PE. Ex:565nm. Em:574nm.

Background

Recommended Dilutions

FC 5 μ l per 10^6 cells in 100 μ l
volume

Immunogen Information

Gene ID

6964

Swiss Prot

P04234

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

TCRD; TRD@; TCRDV1

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Mouse

Isotype

IgG

Purification


Affinity purification

Storage

Store at 2-8°C. Avoid freeze.

Buffer: PBS with 0.09% Sodium azide, 0.2% BSA, pH7.3.

Validation Data



Flow cytometry: 1X10⁶ Human PBMC were surface-stained with ABflo® 450 Rabbit anti-Human/Monkey CD3 mAb (A27177, 5 µl/Test) and PE Mouse IgG1 isotype control (A22878, 5 µl/Test, left) or PE Mouse anti-Human TCR V62 mAb (A27981, 5 µl/Test, right). Cells in the lymphocyte gate were used for analysis.