

ABflo® 500 Mouse IgG 1 isotype control

Catalog No.: A28041

Basic Information

Observed MW

Calculated MW

Category

Primary antibody

Applications

FC

Cross-Reactivity

CloneNo number

AMC50005

Conjugate

ABflo® 500. Ex:410nm. Em:501nm.

Background

Recommended Dilutions

FC ≤0.5 µg per million cells
in 100 µl volume

Immunogen Information

Gene ID

/

Swiss Prot

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Mouse

Isotype

IgG

Purification


Affinity purification

Storage

Store at 2-8°C. Avoid freeze.

Buffer: PBS with 0.09% Sodium azide, 0.2% BSA, pH7.3.

Validation Data



Flow cytometry: 1×10^6 Human PBMC were surface-stained with ABflo® 500 Mouse IgG1 isotype control (A28041, 0.5 μ g, left) or ABflo® 500 Mouse anti-Human CD56 mAb (A26643, 0.5 μ g, right). Cells in the lymphocyte gate were used for analysis.