

Biotin Rabbit anti-Mouse IgA (Fc) mAb

Catalog No.: A28356

Basic Information

Observed MW

Calculated MW

Category

Primary antibody

Applications

FC (intra)

Cross-Reactivity

Mouse

CloneNo number

ARC71353

Conjugate

Biotin

Recommended Dilutions

FC (intra) ≤0.25 µg per million cells
in 100 µl volume

Background

Immunogen Information

Gene ID

/

Swiss Prot

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at 2-8°C. Avoid freeze.

Buffer: PBS with 0.09% Sodium azide, 0.2% BSA, pH7.3.

Validation Data



Flow cytometry: 1×10^6 293T (Transfection) cells were intracellularly-stained with Biotin Rabbit IgG isotype control (A25626, 5 μ l/Test, left) or Biotin Rabbit anti-Mouse IgA (Fc) mAb (A28356, 0.25 μ g, right), followed by APC Streptavidin staining.