

ABflo® 450 anti Mouse CXCR5/CD185 mAb

Catalog No.: A28440

Basic Information

Observed MW

Calculated MW

42 kDa

Category

Primary antibody

Applications

FC

Cross-Reactivity

Mouse

Conjugate

ABflo® 450. Ex:406nm. Em:445nm.

Background

Predicted to be involved in response to X-ray. Located in endoplasmic reticulum.

Recommended Dilutions

FC ≤0.25 µg per million cells
in 100 µl volume

Immunogen Information

Gene ID

12145

Swiss Prot

Q04683

Immunogen

This information is considered to be commercially sensitive.

Synonyms

IKIP

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rat

Isotype

IgG2a, κ

Purification

Affinity purification

Storage

Store at 2-8°C. Avoid freeze.

Buffer: PBS with BSA,pH7.2.