

PE-Cy7 anti-Mouse CD117/c-Kit mAb

Catalog No.: A28479

Basic Information

Observed MW

Calculated MW

109kDa

Category

Primary antibody

Applications

FC

Cross-Reactivity

Mouse

Conjugate

PE-Cy7. Ex:565nm. Em:778nm.

Background

The c-Kit proto-oncogene is the cellular homolog of the transforming gene of a feline retrovirus (v-Kit). The c-kit protein includes characteristics of a protein kinase transmembrane receptor. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

Recommended Dilutions

FC 1 µg per million cells in
100 µl volume

Immunogen Information

Gene ID

16590

Swiss Prot

P05532

Immunogen

This information is considered to be commercially sensitive.

Synonyms

W; Bs; Fdc; Ssm; SCO1; SCO5; SOW3; CD117; c-KIT; Tr-kit; Gsfsc01; Gsfsc05; Gsfscow3

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rat

Isotype

Rat / IgG2b, kappa

Purification

Affinity purification

Storage

Store at 2-8°C. Avoid freeze.

Buffer: PBS with 0.09% Sodium azide and BSA ,pH7.2.