

# Anti-Monkey IgD mAb

Catalog No.: A28758

## Basic Information

---

**Observed MW**

**Calculated MW**

**Category**

Primary antibody

**Applications**

IHC-P

**Cross-Reactivity**

Cynomolgus monkey

**CloneNo number**

ARC5311-01

## Background

---

## Recommended Dilutions

---

IHC-P 1:2000 - 1:8000

## Immunogen Information

---

**Gene ID**

□

**Swiss Prot**

**Immunogen**

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

**Synonyms**

## Contact

---

☎ | 400-999-6126

✉ | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

🌐 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

---

**Source**

Rabbit

**Isotype**

IgG

**Purification**

Affinity purification

**Storage**

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol and 0.05% BSA, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.